

WHISPERGLYDE[®] STATICPOINT™ INSTALLATION STEP BY STEP PROCEDURE FOR EXISTING ROD

WHAT YOU WILL NEED:

STEP LADDER



2 – ADJUSTABLE WRENCHES



TIME REQUIRED: APPROXIMATELY 3-8 MINUTES

DIFFICULTY LEVEL: SIMPLE



SAFETY FIRST! ENSURE THAT YOUR WORK AREA IS CLEAR. FALLING OBJECTS MAY CAUSE INJURIES. ALWAYS HAVE AN ASSISTANT HOLD THE LADDER.

STEP 1

Locate or install a threaded rod adequately supported to the structure. It looks like this:



STEP 2

Confirm that threaded rod matches thread size on your StaticPoint™.



STEP 3

Place a jam nut on the rod before the StaticPoint™. Thread StaticPoint™ onto the threaded rod.



STEP 5

Ensure that there is no more than 2 inches of rod exposed below a lateral brace. Anything longer could result in bending the rod if a sudden heavy load is applied.



STEP 4

Using 2 adjustable wrenches or an end wrench, tighten the jam nut against the StaticPoint™.



STEP 6

The StaticPoint™ “caster” end is intended to provide freedom of circular motion.



If more than 2 inches of rod is unsupported laterally, a “diamond” shaped load ring StaticPoint™ can be purchased to avoid bending the rod.

**TROLLEY MAXIMUM LOAD RATING – 500 LB
DO NOT EXCEED LOAD RATING**



Congratulations on your purchase of the *WhisperGLYDE®* Balance Therapy System. This system was specifically designed and engineered to minimize noise while providing the most efficient trolley movement for the highest quality patient experience and therapy center working environment. Designed to work with *FYZICAL®* balance center *Solostep®* and *Guldmann®* rail systems, the *WhisperGLYDE®* trolley and tether are a direct retrofit. Noise cancelling wheels with soft tyres ensure that the system glides efficiently along the tracks without the harmonic resonance that has been known to elicit complaints from patients and neighboring businesses.

Constructed of aviation grade materials the *WhisperGLYDE®* trolleys also feature low friction bearings and bushings that provide 360 degrees of mobility to patients and durable, trouble free functionality.

WhisperGLYDE® is designed to provide your clients with confidence and peace of mind as they receive their therapy, without the distraction of inconsistent movement of the tether and the high ambient noise that could cause them to lose focus; therapy of mind and body.

We appreciate your support and commitment to *WhisperGLYDE Technology®*. Feel free to contact us at +1.941.360.1150 Extension 1, or info@trixpm.com for sales and technical support.



INSPECTION & MAINTENANCE

- Always read all instructions and warnings contained on the product and packaging before using *WhisperGLYDE®* balance therapy equipment.
- Only use *WhisperGLYDE®* trolleys, tethers and harnesses on approved rail systems: *SoloStep®* or *Guldmann®* for *FYZICAL®*.
- Do not alter or interchange non-matched components.
- Do not disassemble any elements of the *WhisperGLYDE®* trolleys for any reason.
- All balance therapy equipment should be inspected prior to each use.
- This system is not intended to lift the patient; but only to provide support while receiving physical and balance therapy.
- All therapists should be trained in the proper use of balance therapy equipment; including harness application, harness adjustment, tether installation and use.
- Stop using the system if you experience unusual trolley functionality such as intermittent or inconsistent travel and motion. Contact *TRIAx* for technical support.
- Stop using the system if you notice any worn components or looseness in the trolley relative to the rails. Contact *TRIAx* for technical support.
- Follow the recommended installation, inspection and maintenance of the rail system manufacturer.
- Trolleys are "self-braking" under load. They roll freely in normal operation, and do not roll freely when a load is applied.

- Basic care of all safety equipment will prolong the durable life of the unit and will contribute toward the performance of its vital safety function.

HARNESS & TETHER INSPECTION

- Webbing - Look for frayed edges, broken fibers, pulled stitches, cuts, burns, and chemical damage. Do not use if there is any evidence of damaged or compromised webbing straps.
- Connection Points/Loops - Check for distortion, cracks, tears, broken stitching and/or fraying.
- Attachment of Buckles - Inspect for any unusual wear, cracks, frayed or cut fibers, or breaks in the buckle or attachments.
- Carabiners - Inspect to ensure that carabiners are free of distortion in shape and motion. Ensure that secondary locking mechanism is functioning properly and securely tightened.
- Rope Tether - Rotate the tether while inspecting from end-to-end for any fuzzy, worn, broken or cut fibers. Weakened areas from extreme loads will appear as a noticeable change in original diameter. The rope diameter should be uniform throughout, following a short break-in period.

CARE & CLEANING

- Keep rails clean of dirt or other forms of contaminants, and follow all care, maintenance and safety recommend-

- ations of the rail manufacturer.
- Trolleys and related hardware can be cleaned with a mild solution of soap and water. They are self-lubricated and not intended to be lubricated by the user.
- Nylon or Polyester - Remove all surface dirt with a sponge dampened in plain water. Squeeze the sponge dry. Dip the sponge in a mild solution of water and commercial soap or detergent. Work up a lather with a vigorous back and forth motion; then wipe with a clean cloth. Hang freely to dry, but away from excessive heat.
- Drying - Equipment should air dry thoroughly without close exposure to heat, steam, or long periods of sunlight.

WARRANTY & SERVICE

- *WhisperGLYDE®* trolleys have a 3 year warranty against any faults in materials or manufacture. Exclusions from the warranty: normal wear and tear, oxidation, modifications or alterations, incorrect storage, poor maintenance, overload, damage due to accidents, damage due to negligence, and damage from uses for which the product was not designed. *TRIAx* is not responsible for the consequences, direct, indirect or accidental, or any other type of damage befalling or resulting from the use of its products.
- Service replacement products are subject to a refundable core charge.
- For warranty or service, contact *TRIAx* at 941-360-1160 or email info@trixpm.com.